# VICAIR® CENTRE RELIEF

### OFFLOADING PRESSURE IN THE HIGH RISK AREA



### For completely offloading of the rectal area **Adjustable** non-contact area

## Stable positioning

The Vicair Centre Relief is the best choice if you need to completely offload the rectum perineum and coccyx area. By simply removing SmartCells from the centre circle compartment you create a non-contact area. This important feature makes a great difference for persons with (for instance) localized deep pressure ulcers, rectum or tailbone problems. Or after perineal surgery.

The Vicair Centre Relief helps you to sit comfortably and safe, even in the most painful situations.

| Skin protection / Anti decubitus  |             |
|---|-------------|
| Pressure redistribution   | <b>AAAA</b> |
| Shear force reduction   |             |
| Microclimate control  |             |
| Impact & vibration damping  | ***         |
|   |             |
| Positioning   |             |
| Lateral stability   |             |
| Frontal stability   |             |
| Sliding reduction   |             |
| Allows manual adjustments   |             |
| Allows easy transfers   |             |
|   |             |
|   |             |
| Maintenance & Safety  |             |
| Maintenance & Safety  Durability  |             |
|   |             |
| Durability  | ***         |
| Durability<br>Reliability   |             |
| Durability Reliability Functionality when misused   |             |
| Durability Reliability Functionality when misused Low maintenance level   |             |
| Durability Reliability Functionality when misused Low maintenance level   |             |
| Durability Reliability Functionality when misused Low maintenance level Level of hygiene                                      |             |
| Durability Reliability Functionality when misused Low maintenance level Level of hygiene Functionality & Application          |             |
| Durability Reliability Functionality when misused Low maintenance level Level of hygiene  Functionality & Application Comfort |             |

Low seat to floor height Auto adjustments

### **USED WITH**

Rectum problems, tailbone problems, localized deep pressure ulcers, after perineal surgery or difficult childbirth

| TECHNICAL SPECIFICATIONS |           |  |
|--------------------------|-----------|--|
| Height                   | 10 cm     |  |
| Average weight           | 750 grams |  |
| Max. user weight         | 250 kg    |  |
|                          |           |  |

#### **CUSHION SIZES**

View the Morton & Perry comparison chart for available sizes and product codes.

070069-ACR-EN-REV04/2019