VICAIR® ACTIVE 02

FOR THE ACTIVE WHEELCHAIR USER











PRESSURE ULCER TREATMENT

Perfect for users who transfer frequently

Excellent lateral and frontal stability

Features a nifty storage pouch for personal items

The Active O2 offers the unique combination of ergonomically shaped foam and Vicair SmartCells. The four rear SmartCell compartments provide a high level of comfort and skin protection, and reduce pressure on the sitting bones and tailbone area. The foam front makes transfers easy and feels steady and improves wheelchair handling. The foam inlay is protected by an easy-toclean fluid-resistant liner. The rest of the cushion is 100% machine washable and breathable for optimal hygiene and perfect temperature & moisture regulation. The Active O2 has an easy accessible storage pouch for precious belongings (keys, phone, wallet). In short: everything you need to maintain your active lifestyle!

PRODUCT PERFORMANCE	
Skin protection / Anti decubitus	***
Pressure redistribution	***
Shear force reduction	***
Microclimate control	
Impact & vibration damping	***
Positioning	***
Lateral stability	***
Frontal stability	***
Sliding reduction	***
Allows manual adjustments	***
Allows easy transfers	***

Maintenance & Safety	***
Durability	***
Reliability	4444
Functionality when misused	
Low maintenance level	AAAA
Level of hygiene	AAAA

Functionality & Application	***
Comfort	AAAA
Low cushion weight	
Available size range	AAAA
Maximum user weight	AAAA
Low seat to floor height	AAAA
Auto adjustments	

USED WITH

Spinal cord injuries, active users with a high transfer frequency

TECHNICAL SPECIFICATIONS		
6 and 9 cm		
1250 grams		
200 kg		
60°C (without foam)		

CUSHION SIZES

View the Morton & Perry comparison chart for available sizes and product codes.

070069-OCT-EN-REV04/2019