

**ULTRA-LIGHTWEIGHT, EASY TO USE AND COMFORTABLE**

Compartment barrier  
gently limits sliding forward

Equal pressure distribution  
on the seating surface

Auto-adjusting valve



**Air valve automatically adapts the cushion to your body shape and weight**

**Two compartment layout limits sliding**

**Easy to clean and disinfect**

The Vicair Liberty is an easy to use and comfortable (wheelchair) cushion. It is designed for users who require comfort in combination with a low to medium level of support and pressure ulcer prevention. The cushion is ready-for-use, rain proof and with two compartments it provides an even pressure distribution. This makes it an ideal cushion for occasional use in a large variety of situations.

This Vicair cushion is extremely lightweight.

**PRODUCT PERFORMANCE**

|                                  |       |
|----------------------------------|-------|
| Skin protection / Anti decubitus | ▲▲▲▲△ |
| Pressure redistribution          | ▲▲▲▲△ |
| Shear force reduction            | ▲▲▲▲△ |
| Microclimate control             | ▲▲▲▲△ |
| Impact & vibration damping       | ▲▲▲▲△ |

|                           |       |
|---------------------------|-------|
| <b>Positioning</b>        | ▲▲▲▲△ |
| Lateral stability         | ▲▲▲▲△ |
| Frontal stability         | ▲▲▲▲△ |
| Sliding reduction         | ▲▲▲▲△ |
| Allows manual adjustments | ▲▲▲▲△ |
| Allows easy transfers     | ▲▲▲▲△ |

|                                 |       |
|---------------------------------|-------|
| <b>Maintenance &amp; Safety</b> | ▲▲▲▲△ |
| Durability                      | ▲▲▲▲△ |
| Reliability                     | ▲▲▲▲△ |
| Functionality when misused      | ▲▲▲▲△ |
| Low maintenance level           | ▲▲▲▲△ |
| Level of hygiene                | ▲▲▲▲△ |

|  |       |
|--|-------|
| <b>Functionality &amp; Application</b> | ▲▲▲▲△ |
| Comfort                                | ▲▲▲▲△ |
| Low cushion weight                     | ▲▲▲▲△ |
| Available size range                   | ▲▲▲▲△ |
| Maximum user weight                    | ▲▲▲▲△ |
| Low seat to floor height               | ▲▲▲▲△ |
| Auto adjustments                       | ▲▲▲▲△ |

**USED WITH**

Osteoarthritis/arthritis, rheumatism, spondylitis, pain relief, mobility problems, in a geriatric environment.

**TECHNICAL SPECIFICATIONS**

|                  |           |
|------------------|-----------|
| Height           | 5 cm      |
| Average weight   | 350 grams |
| Max. user weight | 150 kg    |

**CUSHION SIZES**

View the Morton & Perry comparison chart for available sizes and product codes.