



**CLICK ON ANY IMAGE
TO GO TO THE PRODUCT PAGE**



VICAIR

NOVIS

NOVIS

NOVIS

NOVIS

NOVIS

Maxxcare

MODEL	Aircell	Better Living	ProCair	ProCair Prime	ProCair Plus	CairMax Duo	Heel Pro Evolution
Product Type	Aircells in compartments	Alternating Air Mattress	Alternating Air Mattress	Alternating Air Mattress	Alternating Air Mattress	Hybrid Aircells with Foam	Aircells in compartments
INDICATION FOR USE							
Ideal for	Ideal for those who have sleep, pain, comfort and positioning issues. Pressure redistribution is achieved through immersion. A non-powered static system no daily maintenance required.	Mattress replacement for home and residential pressure care. Easy and quick setup. Automatic pressure setting and cell in cell to enhance user safety and comfort. Very economically priced.	For all healthcare environments. With automatic weight sensing pressure care. Very easy setup. Economically priced.	A mattress with automatic weight sensing and an enclosed, out of sight pump. A handcontrol allows the user to operate all therapy functions while in bed. With a foam underlay that can be positioned under or on top of the cells as required.	Homecare and professional mattress with the following safety benefits: Auto sensing, fowler boost, bottoming out protection, heel zone, exit alarm and battery backup.	Ideal for those who require improved comfort and potentially changing needs. Hybrid mattress, with step up (with pump) and down (without pump) therapy. Very economical priced.	For those at risk of a heel injury or with an existing pressure injury. Superior technology and benefits versus foam and fibre products.
PRESSURE INJURY RISK - Indication							
Low- Very High	●	●	●	●	●	●	●
SPECIFICATIONS							
SWL Kg	<250Kg	40-180Kg	40-200Kg	30-200Kg	40-220Kg	40-230Kg	N/A
CPR	Yes, without the use of a board	Quick Twist	Quick Twist	Quick Twist	Quick Twist	No	N/A
Pump	N/A	Yes	Yes	Yes	Yes	Yes, optional	N/A
Hybrid conversion option	N/A	N/A	N/A	Yes	N/A	Yes	N/A
Auto weight sensing for pressure	Aircells, self adjust	Yes	Yes	Yes	Yes	Yes	N/A
Fowler Boost for backrest raise	N/A	No	Manual	No	Automatic	No	N/A
Bottoming out protection	Yes	Yes	Yes	Yes	Yes	Yes	N/A
Warranty	3 Year	1 Year	3 Year	3 Year	3 Year	2 Year	6 Months
AVAILABLE SIZES Width x Length x Depth - product order codes							
Single 88x200x20cm	XVARM415	XNAPMBL-R01	XNAPMPC-R01	XNAPMC-R03	XNAPMPC-R02	XNFAMCM-R01	
King Single 1050x200x20cm	XVARMKS415		XNAPMPC-R01K		XNAPMPC-R02K	XNFAMCM-R01K	
Custom on application	Yes						
Small							XVHEELSML
Standard							XVHEEL
Extra Large							XVHEELXL

Maxxcare
Measuring sizing for:
Heel Pro Evolution



Small: < 38cm / 15"
Standard: 38-43cm / 15-17"
XL: ≥ 43cm / 17"



CLICK ON ANY IMAGE
TO GO TO THE PRODUCT PAGE

VICAIR®

Mattress 415

Static aircells in multiple compartments. Ideal for sleep, pain, comfort and positioning.

Versus other static air systems:

- > A Mattress replacement not an overlay.
- > Lower maintenance and improved compliance. No daily monitoring, inflation or pumps required.
- > Individual compartments are adjustable in cell density if required.



Max SWL 350Kg

Also available in wider 1050mm King Single and custom sizes.

novis | BetterLiving™

Alternating air mattress replacement

Intended for homecare and residential care.

- > Very economically priced
- > Automatic weight and pressure setting
- > Cell in cell for enhanced user comfort and safety (prevents bottoming out)
- > Easy to operate and maintain



novis | CairMax Duo

Hybrid foam air mattress replacement

Ideal for those who require more comfort.

Duo equals two separate support options in one mattress:

- Option 1** Without the pump - provides surface pressure distribution. (foam only)
- Option 2** Attach alternating pump for dynamic therapeutic pressure relief (foam + air).

Features

- > Automatic pressure and weight setting
- > Easy to operate and maintain.




novis
ProCair range
Compare the benefits



	ProCair	ProCair Prime	ProCair Plus
AutoCair Automatic weight sensing technology that delivers fail safe patient comfort levels. No more concerns over correct pressure level settings. Manual adjustments can still be made.	✓ Cyclic	✓ Cyclic	✓ Real-time
Quicker Release CPR access Easier quick twist mechanism located at the head of the mattress. Makes it easier to rapidly deflate the mattress for emergency CPR requirements.	✓ Quick Twist	✓ Quick Twist	✓ Quick Twist
Fowler Boost Adjusts mattress pressure to provide increased support under the sacral area when backrest is raised.	✓ Manual		✓ Automatic
SafetyCair The ultimate protection from 'bottoming out', alternating air cells contain a lower 'safety cell' which remains fully inflated to provide greater surface stability and patient safety.	✓ Air Cell	✓ Foam Insert	✓ Air Cell
HeelCair Independent heel zone made up of narrower cells. Designed for individualised therapy to the 'at risk' heel and lower leg area.			✓ Disconnectable
ClimateCair ClimateCair cells are constructed from thermoplastic polyurethane (TPU) material. Not only softer to the touch, the cells offer enhanced heat retention qualities, to maximise patient comfort.			✓
CairAlert Be alerted when your patient leaves the bed with a warning alarm that, when activated, will signal the carer that the patient has left the mattress. Ideal for patients with cognitive dysfunction.			✓
FailSafe An inbuilt Back Up Battery unit provides emergency back up power for continuous operation up to 5 hours. Also ideal for alternating therapy during extended patient transportation periods.			✓
King Single Wider surface for therapy and positioning.	✓		✓

Maxxcare®
Heel Pro Evolution

Prevents and treats heel pressure injuries



2nd Most common pressure injuries are on the heel