

Sealed base

A TPU laminated nylon replaces the unlaminated, permeable base allowing for wipe down of surfaces and greater infection control.



Larger CPR valve

Faster deflation in the event of an emergency and faster packing for when the mattress is stored.



Exposed umbilical

Easy wipe down, increased infection control without nylon sleeve.

a. The sleeve did not provide an impermeable seal at the umbilical connector. Clinical recommendations from healthcare professionals indicate that umbilical cords are now easily cleaned and maintained.

b. The umbilical cords are no longer welded together. Welded joins are an infection control risk. There are now 2 independent cords connecting to the pump for optimal cleaning and disinfecting.

