

Independence Chairs



CLICK ON ANY IMAGE
TO GO TO THE PRODUCT PAGE



Morton & Perry
Healthcare Equipment Solutions



MODEL	VELA Tango 50	VELA Tango 100	VELA Tango 200EI	VELA Tango 500	VELA Tango 500EI	VELA Tango 310E
Item Code	XVE6010	XVE6021	XVE6149	XVE6045	XVE6155	XVE5106
Standard	All chairs are stable, with a central brake, height adjustable (either by gas or electric). Seat depth is increased by backrest positioning. Backrest height adjustable. Armrests adjustable up/down and for width.					
Additional Features	> Low backrest supports lower back and allows free movement of the back, shoulders and arms > Includes a footrest as standard	> 360° rotation seat > Low backrest supports lower back and allows free movement of the back, shoulders and arms	> Higher Backrest > Larger Seat W47 x D49cm > Electric height adjustment	Same as Tango 200 with addition of > 2 Part ALB Backrest > Adjustable back rest lumbar support > Included seat rotation 90° LH and RH > Fold down armrests	Same as Tango 500 with addition of > Electric height adjustment	> Robust chair for 300kg user weight > Electric height adjustment to assist user to get up from the chair > Includes 4XL seat (W: 72xD: 50 cm) and a CX backrest (W: 38xH: 43 cm)
SPECIFICATIONS						
Seat lift	15cm - Gas Lift	15cm - Gas Lift	20cm - Electric Lift	15cm - Gas Lift	20cm - Electric Lift	15cm - Electric Lift
Seat height - Fixed Position 1*	50 - 65cm	50 - 65cm	50 - 65cm	50 - 65cm	50 - 65cm	42 - 57cm
Seat height - Fixed Position 2*	48 - 63cm	48 - 63cm	48 - 68cm	48 - 63cm	48 - 68cm	
Seat - width x depth	44 x 44cm or large: 47 x 49cm	44 x 44cm	47 x 49cm (18.5 - 19.3 inches)	44 x 44cm (17.3 - 17.3 inches)	44 x 44cm (17.3 - 17.3 inches)	72 x 50cm (28.3 - 19.7 inches)
Seat tilt	-	6° Forwards, 3° backwards	6° Forwards, 3° backwards	90° Forwards, 22° backwards	10° Forwards, 7° backwards	-
Seat rotation	-	Standard 360°	-	90° to the right / 90° to the left	90° to the right / 90° to the left	-
Sitting depth	37 - 49cm (14.6 - 19.3 inches)	37 - 44cm (14.6 - 17.3 inches)	33 - 49cm (13 - 19.3 inches)	42 - 48cm (16.5 - 18.9 inches)	42 - 48cm (16.5 - 18.9 inches)	46 - 56cm (18 - 22 inches)
Backrest tilt	10° Forwards, 13° backwards	10° Forwards, 13° backwards	20° Forwards, 6° backwards	4° Forwards, 22° backwards	4° Forwards, 22° backwards	-
Height adjustment of backrest	10cm (3.9 inches)	10cm (3.9 inches)	13cm (5.1 inches)	4cm (1.6 inches)	4cm (1.6 inches)	0 - 15cm
Armrest height	Adjustable 0 - 24cm (0 - 9.4 inches)	Adjustable 0 - 24cm (0 - 9.4 inches)	Adjustable 0 - 24cm (0 - 9.4 inches)	Adjustable & foldable 0 - 22cm (0 - 8.7 inches)	Adjustable & foldable 0 - 22cm (0 - 8.7 inches)	Adjustable 0 - 24cm (0 - 9.4 inches)
Distance between armrests	Adjustable 37 - 51cm (14.6 - 20 inches)	Adjustable 37 - 51cm (14.6 - 20 inches)	Adjustable 41 - 55cm (16.1 - 21.7 inches)	Adjustable 42 - 54cm (16.5 - 21.2 inches)	Adjustable 42 - 54cm (16.5 - 21.2 inches)	Adjustable 70 - 75cm (27.5 - 29.5 inches)
Foot rest - Fold up	Yes, standard	Optional	Optional	Optional	Optional	-
Chair's weight	20kg	20kg	29kg	26kg	33kg	42kg
Max. user weight	125kg	160kg	160kg	160kg	160kg	300kg

* Seat Height depends on the fixed height setting of the chair.
There are two fixed settings.

Further modifications and accessories are available on request





Kitchen



Dining



Ironing



Work space

Greater independence and dignity

Ideal for the following conditions:

- > Spinal disorder
- > Mobility impairment
- > Walking impairment
- > Arthritis
- > Multiple sclerosis
- > Cerebral palsy
- > Neurology
- > Muscular dystrophy
- > Amputation
- > Parkinsonism
- > Back pain
- > Stroke
- > Restricted growth

Optional Accessories



Strolling Bracket

Improves the working environment of the care provider, e.g. when the user is moved sitting in the chair. Recommended to use with foot rest.



Foot rest - Fold up

The feet can rest while sitting or moving around. Can also be used as a step when getting in and out of the chair.



Larger CY Seat

The CY seat (W:47xD:49 cm) is a comfortable, padded seat with a rounded front. The user enjoys good, supportive sitting comfort, and the seat design also ensures good blood flow to the legs.



T-grip for break handle

People with pains and reduced joint mobility may benefit from the T-grip. It reduces twists in the body and back, when using the brake on the chair. To increase the height, a brake extender is also available.

Tango a Unique Chair

More functional than a wheelchair or walker for indoor tasks

Support & Comfort

Backrests - All Tango models are available with a varied selection of backrests so that the chair can be adapted to the individual user. Each back can tilt and is height adjustable.

Ergonomic

Seat options - The adjustment of the seat height, depth and tilt plus adjustable backrests and armrests ensures the optimum ergonomic sitting position.

Height Adjustable - All chair models are available with manual or electric height adjustment. The electric height adjustment can be operated from the left or right hand.

Easier Transfers

Armrests - The soft armrests are height and width adjustable. They can also lower down to seat height for easier transfers to seat, bed, toilet or wheelchair.

Safe

Central Brake
The Vela chair has a unique braking system that ensures the user a stable and comfortable workplace.

Walkable

Directional Wheels - The chair is equipped with four smoothly running directionally stable wheels that allow the chair to be "walked" around so that the user can get around at home unassisted.

Footrest

Foldable Footrest - The feet can rest while sitting or while being moved around.



Award Winning Design

VELA Tango 500 with ALB backrest (pictured above) won the prize for "Best New Design" at the Scandinavia Health & Rehab 2016 fair in Copenhagen.



Increased Reach



Personal Hygiene



Transfers



Cooking